



Piezo devices



Piezo devices

Using the 2 devices



Connect the 2 devices
through
Y splitter
To iPhone/iPad




Y splitter cable is specific for piezo devices
Never use it with Plasma balls
Never use the Plasma Ball splitter for piezo devices




Piezo devices

Using a single device






Connect the device
Directly
To iPhone/iPad

Piezo devices




Pink side should be placed to your skin

VOLUME of iPhone/iPad (small buttons on the left/right of your device)
Should be placed at **50%**




COIL & PLASMA

Use *COIL & PLASMA* option



Piezo devices

Using the 2 devices



To treat a specific part of the body


Seems to boost stem cells

Can reduce pain

Piezo devices

Using the 2 devices


For example:
You have pain
or need to stimulate your liver



The screenshot shows a software interface with a list of frequencies and their corresponding organs. A pink arrow points to the list.


Frequency	Organ
60	Colon - Asc & Desc; Sacrum; Uterus; Ovaries; Small Intestine;
61	Uterus; Ovaries; Small Intestine; Transverse Colon; Stomach; Belly;
62	Stomach; Pancreas; L5=>L1;
63	Spleen; Gall Bladder; L5=>L1;
64	Spleen; Gall Bladder; Liver;
65	Liver; Kidneys; Adrenal Glands; T12=>T9;
66	Chest; Breasts; T12=>T9;
67	Chest; Breasts; Diaphragm;
68	Diaphragm; Aorta - Descending;

Other interface elements include: PLAY, STOP, BEAT 0.1 Hz, SHIFT 180°, 10 SEC, Frequencies TABLES, (C 1) Freq Selected: 169.3 Hz, Freq. Casc. 1, LOAD Seq., NEW Seq., COIL-PLASMA-THERAPI SINE - Phi snds CASCADE UP, 5 & => 10 Spleen; Gall Bladder; Liver;, CASCADE of SINGLE Freq., PAIR of RIFE Freq., and a list of frequencies: 1 Hz-Toe, 59 Hz-Tendon, 62 Hz-Tibia, 121 Hz-Anus Coccyx, 130 Hz-Prostate, 254 Hz-Lungs, 264 Hz-Lungs, 513 Hz-Cerebellum, 554 Hz-Pineal Gl., 896 Hz-Crown Chk.




The screenshot shows a software interface with a diagram of a human body. A pink arrow points to the diagram.

Other interface elements include: PLAY, STOP, BEAT 0.1 Hz, SHIFT 180°, 10 SEC, Frequencies TABLES, (C 1) Freq Selected: 169.3 Hz, Freq. Casc. 1, COIL-PLASMA-THERAPI SINE - Phi snds CASCADE UP, 5 & => 10 Spleen; Gall Bladder; Liver;, COIL & PLASMA, and a list of frequencies: 1 Hz-Toe, 59 Hz-Tendon, 62 Hz-Tibia, 121 Hz-Anus Coccyx, 130 Hz-Prostate, 254 Hz-Lungs, 264 Hz-Lungs, 513 Hz-Cerebellum, 554 Hz-Pineal Gl., 896 Hz-Crown Chk.



Place one device
on your belly (on your liver)
The other one on your back




The 2 devices are in front of each other






Select ***ONLY SHIFT***

Than

Press ***PLAY***
For 1-2 minutes
(Not more than 2 minutes)






Than, select ***BEAT 0.1Hz***

And

Press ***PLAY***
For a few minutes

Piezo devices

Using a single device




Connect the device
Directly
To iPhone/iPad

Place the device on your lower back part (sacrum - Hara)

You will need to create and play the frequencies of your chakras

How to create a sequence of body frequencies (1)



To create a 10 minutes sequence, select **10 SEC**
To create a 20 minutes sequence, select **20 SEC**

How to create a sequence of body frequencies (2)

Select the first frequency

Nbr
in the sequence

Is
Selected

64 - Spleen; Gall Bladder; Liver;	165.33 Hz
65 - Liver; Kidneys; Adrenal Glands; T12=>T9;	176.35 Hz
[01-10s-A0.3-SiUP5] 66 - Chest; Breasts; T12=>T9;	185.01 Hz
67 - Chest; Breasts; Diaphragm;	196.23 Hz
68 - Diaphragm; Aorta - Descending;	211.67 Hz

COIL-PLASMA-THERAPHI
SINE - Phi snds CASCADE
UP
5 & => 10
Prostate; Appendix; Bladder; Colon -
Asc & Desc;

	[01-10s-A0.3-SiUP5] 66 - Chest; Breasts; T12=>T9;	185.01 Hz
--	---	-----------

10
SEC

The BEAT time can be different for each frequency of the sequence
If so, the BEAT time should be changed
before selecting a frequency in the sequence

How to create a sequence of body frequencies (3)

Select the frequencies needed in the sequence

64 - Spleen; Gall Bladder; Liver;	165.55 Hz
65 - Liver; Kidneys; Adrenal Glands; T12=>T9;	176.35 Hz
[01-10s-A0.3-SiUP5] 66 - Chest; Breasts; T12=>T9;	185.01 Hz
[02-10s-A0.3-SiUP5] 67 - Chest; Breasts; Diaphragm;	196.23 Hz
[03-10s-A0.3-SiUP5] 68 - Diaphragm; Aorta - Descending;	211.67 Hz
[04-10s-A0.3-SiUP5] 69 - Aorta - Descending; T8=>T6;	228.33 Hz
[05-10s-A0.3-SiUP5] 70 - Heart; T8=>T6;	242.18 Hz
[06-10s-A0.3-SiUP5] 71 - Heart; Lungs;	254.01 Hz
[07-10s-A0.3-SiUP5] 72 - Lungs;	264.59 Hz
73 - Lungs; Thymus Gland;	273.99 Hz

It is possible to deselect and go back in the sequence.

Not possible to deselect inside the sequence

Then
SAVE the Seq



Enter the
name
of the Seq

Record Name
Insert the Name,
related to this record.
(Max 10 characters)
If empty =

OK

You can also
leave
a comment

Annotate record
Leave a note to identify your record,
if you want to find it later.

OK

For Chakras sequence

Select from the list

Fréquences number:

05

14

41

59

68

77

86


95

96



51 - Knee;	87.21 Hz
52 - Knee - Ligaments of;	92.50 Hz
53 - Knee - Kneecap; Back of Knee;	97.02 Hz
54 - Thigh;	101.06 Hz
55 - Thigh;	104.65 Hz
56 - Thigh - Top of; Testicles; Perineum;	108.99 Hz
57 - Perineum; Testicles; Anus;	114.34 Hz
58 - Anus; Coccyx; Rectum; Sigmoid Colon; Prostate;	121.28 Hz
59 - Prostate; Appendix; Bladder; Colon - Asc & Desc;	130.82 Hz

How to PLAY a sequence of body frequencies (1)



Select the Seq
in the table

How to PLAY a sequence of body frequencies (2)

File Selected
heart 2020-01-05

NOTE:
seq up

Nbr in Seq:7

66, 67, 68, 69, 70, 71, 72

[BACK to the LIST](#)

[LOAD Sequence](#)

PLAY **STOP** **BEAT** Turning & Stop **10 SEC** **BODY** Frequency

AMPLI 0.3 EDIT Seq

Freq Selected: 130.82 Hz

COIL-PLASMA-THERAPHI
SINE - Phi snds CASCADE
UP

5 & => 10
Prostate; Appendix; Bladder; Colon -
Asc & Desc;

heart 2020-0
Sequence Nbrs :
66, 67, 68, 69, 70,
71, 72

SEQUENCE ON

Prostate;
Appendix;
Bladder;
Colon - Asc & Desc;



PLAY **STOP** **BEAT** 0.25 Hz **10 SEC** **BODY** Frequency

AMPLI 0.3 EDIT Seq

Freq Selected: 185.01 Hz

COIL-PLASMA-THERAPHI
SINE - Phi snds CASCADE
UP

5 & => 10
Chest; Breasts; T12=>T9;

heart 2020-0
Sequence Nbrs :
66, 67, 68, 69, 70,
71, 72

SEQUENCE ON
BEAT Turning AND STOP=>ON
Seq Nbr [66]

Chest;
Breasts;
T12=>T9;




When playing a sequence, the BEAT is automatically selected to
NOT possible to select a single BEAT frequency

BEAT
Turning & Stop

Piezo devices

Using a single device




Place the device on your lower back part (sacrum - Hara)
And play the frequencies sequence you have defined for your chakras

It is very important to be **present** (your mind connected)
to the session you are receiving. **Get rid of all other (bad) thinking.**

Quite sensitive persons can feel through the session
that something is happening in the different parts (chakras) of the body:
it can be tensions appearing/disappearing; some vibrations; or any kind of feelings,...

Piezo devices

Using a single device



The device placed on your lower back part (sacrum - Hara)
And playing the frequencies sequence you have defined for your chakras

It is very important to be **present** (your mind connected)
to the session you are receiving. **Get rid of all other (bad) thinking.**

Do not worry if you do not feel. It can need a lot of practice.

In this case, practice
Microcosmic orbit breath during the session

https://en.wikipedia.org/wiki/Microcosmic_orbit

